FOAMING TEST APPARATUS ASTM D 892 - IP 146 - ISO 6247

The apparatus consists of a frame supporting one or two Pyrex jars and a control box with all the electronics. Thanks to the use of mass flow controllers, it's now possible to easily and automatically measure out the correct amount of air to be blown into the samples. Models AD0892-500 and -510 have a series of knobs and displays to adjust the flow while model AD0892-520 is equipped with a Panel PC with touch screen interface that permits also sequencing the various phases of the test with buzzers alerting the user and screen messages to inform the user about the action required.



AD0892-510

- Powder-coated steel frame supporting one (AD0892-500) or two Pyrex jars (models -510 and 520).
- Anodized aluminum jar covers, each with two holes for the graduate test cylinders: locking clamps to hold the cylinders in position and thermometer jacket for ASTM 12C thermometer.
- One set of rubber stoppers with air tubes (inlet and outlet).
- Electric bath stirrer with shaft and propeller, stainless steel.
- Stainless steel heaters.
- Coils for both air circulation and connection to a chiller in case tap water temperature exceeds 20°C.
- Electronic mass flow controllers with 0-115 ml/min scale and 0.1% F.S. precision, two in AD0892-500 and four in -510 and -520.
- Model AD0892-500 and -510 come with a traditional aluminium control panel with switches to start flow and a knob plus a display to adjust it: a timer permits to stop blowing automatically after 5 min (selectable).
- Model AD0892-520 is equipped with an 8.4" Panel PC that allows sequencing
 of the entire test procedure: messages on the screen and buzzers guide the
 operator through the analisys. Calibration, diagnostic and data storage
 routines.
- Temperature measuring elements: PT100 RTD. Operating range: from ambient to 100°C with +/- 0.1°C accuracy.
- One Ethernet, two USB and one serial port on model AD0892-520.
- Safety transparent shield surrounding 93.5°C bath.
- Safety devices against overheating and low level.
- Easy access control box that can be removed in case of malfunctionment for easy transportation.
- English written user's guide and installation instructions.
- For 220V / 50Hz connections. Power consumption 1800 W (AD0892-500 and 2500 W (-510 and -520)
- Dimensions (I x w x h): mm 550 x 460 x 720 mm, 30 kg (AD-0892-500) mm 1050 x 460 x 720, kg 55 (AD0892-510 and -520)
- CE marked.

AD0892-500 Apparatus, single-unit AD0892-510 Apparatus, twin-unit

AD0892-520 Apparatus, twin-unit with touch-screen interface

ACCESSORIES

AD0892-A20 Electronic air volume counter (replaces wet test meter)

AD0892-A21 Certificate for the volume counter

AD0892-A10 Air pump

AD0892-A03 Glass drying tower CAL001 PT100 simulator

CAL003 Official Certificate for Pt100 simulator

CONSUMABLES

AD0892-C00 Graduate pyrex glass cylinder
AD0892-C01 Rubber stopper with air inlet and outlet
AD0892-C02 Norton certified spherical Alundum gas diffuser stone

AD0892-C03 Mott certified stainless steel cylindrical gas diffuser TA012C-N00 ASTM 12C thermometer (-20°C/+102°C)

THO 120 TWO TOTAL 120 MICHIGATION CLOT (20 MICHIGAN CLOT)

AD0892-520

Specifications may vary without notice.

The apparatus is supplied bare without glassware (if any), accessories and consumables.

